

Microscope Kit

WHAT'S IN THE KIT

- 2 microscopes
- 2 sets of 48 prepared slides: animals, bugs, materials, plants



QUICK-START GUIDE

1. Start at lowest magnification and center the slide.
2. Focus with the coarse knob, then the fine knob.
3. Sketch what you see — drawing forces careful observation.
4. Let students choose slides; curiosity drives the station.

FULL LESSONS IN THIS GUIDE

GRADES K-2	Tiny Worlds Drawing Lab	page 2
GRADES 3-5	Mystery Slide Detective	page 3
GRADES 6-8	Compare & Classify	page 4

PRO TIP: Two scopes work best as a rotation station with a sketch-journal task for waiting students.

Microscope Kit

LESSON · GRADES K-2

Tiny Worlds Drawing Lab

LEARNING OUTCOMES

- Students will use a magnification tool to observe details
- Students will record observations through drawing
- Students will compare what tools reveal versus unaided eyes

MATERIALS

- Microscopes (2 stations) + prepared slide sets
- Two-panel drawing sheet: 'my eyes' vs. 'the microscope'
- Crayons/pencils
- Station rotation timer

INSTRUCTIONS

1. Model slide handling and the focus knob gently — teacher sets magnification.
2. Students first draw the slide subject as their eyes see it (tiny!).
3. At the microscope, draw the same subject as the lens reveals it.
4. Rotate: half the class draws from observation books while half uses scopes.
5. Gallery: post the two-panel drawings side by side.
6. Circle: what other tools help people see, hear, or reach beyond their bodies?

ACTIVITY

Tiny Worlds Drawing Lab — the same bug drawn by eye and by lens, side by side. Tools-extend-our-senses is the whole lesson, and their own drawings prove it.

Microscope Kit

LESSON · GRADES 3-5

Mystery Slide Detective

LEARNING OUTCOMES

- Students will observe systematically and record evidence
- Students will use clues to identify unknown specimens
- Students will defend a conclusion with observed details

MATERIALS

- Microscopes + slide sets with labels covered
- Detective sheet: sketch, features, hypothesis
- Clue cards per specimen (habitat, diet, category hints)
- Reveal envelopes

INSTRUCTIONS

1. Cover slide labels in advance; number each mystery slide.
2. Detectives observe, sketch, and list three concrete features per slide.
3. Clue cards unlock after the sketch is done — evidence before hints.
4. Teams commit to identifications with a written justification.
5. The reveal: open envelopes, compare, and honor near-misses with good reasoning.
6. Debrief: which features were most diagnostic? Which misled?

ACTIVITY

Mystery Slide Detective — observe first, hypothesize second, reveal last. Requiring sketches before clue cards makes evidence-based reasoning the only path to the answer.

Microscope Kit

LESSON · GRADES 6-8

Compare & Classify

LEARNING OUTCOMES

- Students will document cell-level differences between specimen types
- Students will build a classification argument from observations
- Students will use comparative evidence in scientific writing

MATERIALS

- Microscopes + plant and animal slide selections
- Comparison documentation sheet (labeled sketches per specimen)
- Classification argument frame: claim, evidence, reasoning
- Reference images for verification

INSTRUCTIONS

1. Teams observe paired specimens: plant tissue vs. animal tissue slides.
2. Documentation standard: labeled sketch + magnification noted + features listed.
3. Build a comparison table across all observed specimens.
4. Write the classification argument: how do you tell plant from animal at this scale?
5. Peer review: does each claim trace to a documented observation?
6. Test the rule: apply it to a new unlabeled slide.

ACTIVITY

Compare & Classify — from labeled sketches to a defensible classification rule, tested against an unknown slide. Foundational life-science practice with their own evidence chain.